

Anti-Mouse Cd59A Monoclonal Antibody () For Flow Cytometry - O55186

Item Number: CM0009487

Introduction

Rabbit monoclonal antibody against mouse CD59a, conjugated to . Validated for flow cytometry. CD59a inhibits complement membrane attack complex formation, protecting cells from lysis. Cross-reacts with mouse. Protein weight ~14kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD59a mAb
Remarks/Alias	Cd59; MACIF; MAC-IP
Species	Mouse
Gene ID	12509 (mouse)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-96 of mouse CD59a (NP_031678.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	O55186
Protein Weight	14kDa
Shipping	Ice bag