

Anti-Mouse Cd69 Monoclonal Antibody () For Flow Cytometry - P37217

Item Number: CM0010719

Introduction

-conjugated rabbit monoclonal antibody specific for mouse CD69, an early activation antigen on T cells, B cells, NK cells. Validated for flow cytometry (FC), ideal for immunology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD69 mAb
Remarks/Alias	AIM; VEA; 5830438K24Rik
Species	Mouse
GeneID (Human)	969
GeneID	12515
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 64-199 of mouse CD69 (NP_001028294.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P37217
Protein Weight	22kDa
Shipping	Ice bag