

Anti-Mouse Cd79B Monoclonal Antibody () For Flow Cytometry - P15530

Item Number: CM0010858



Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD79b (Ig-beta) for flow cytometry. Ideal for B-cell research, BCR signaling studies, and immunophenotyping. Specific to mouse CD79b extracellular domain.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD79b mAb
Remarks/Alias	B29; Igb; Igbeta; Ig-beta
Species	Mouse
GeneID (Human)	974
GeneID	15985
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-158 of mouse CD79b (NP_032365.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P15530
Protein Weight	25kDa
Shipping	Ice bag