

Anti-Mouse Cd79B Monoclonal Antibody For Flow Cytometry - P15530

Item Number: CM0008537



Introduction

-conjugated rabbit anti-mouse CD79b monoclonal antibody for flow cytometry. Targets B-cell receptor complex beta chain (Ig-beta, B29), validated in mouse models. Ideal for IVD research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD79b mAb
Remarks/Alias	B29; Igb; Igbeta; Ig-beta
Species	Mouse
GeneID (human)	974
GeneID	15985
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 26-158 of mouse CD79b (NP_032365.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P15530
Protein Weight	25kDa
Shipping	Ice bag