

Anti-Mouse Cd80 Monoclonal Antibody () For Fc - Q00609

Item Number: CM0008093



Introduction

-conjugated anti-mouse CD80 monoclonal antibody for flow cytometry. Targets the B7-1 costimulatory molecule involved in T-cell activation and immune regulation. Suitable for immunology research.

[Learn More](#)

Parameter	Value
Product Name	anti-Mouse CD80 mAb
Remarks/Alias	B71; Ly53; TSA1; Cd28l; Ly-53; MIC17
Species	Mouse
Gene ID (Human)	941
Gene ID	12519
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 38-245 of Mouse CD80.
Source	Armenian Hamster
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q00609
Protein Weight	34kDa
Shipping	Ice bag