

Anti-Mouse Cd84 Monoclonal Antibody () For Flow Cytometry - Q18Pi6

Item Number: CM0005432

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD84 (SLAMF5, CDw84). Designed for flow cytometry applications. Recognizes mouse CD84 with a predicted molecular weight of ~37 kDa. SWISS: Q18PI6.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD84 mAb
Remarks/Alias	CDw84, SLAMF5, A130013D22Rik
Species	Mouse
GeneID (Human)	8832
GeneID (Mouse)	12523
Immunogen	Recombinant fusion protein containing amino acids 22-221 of mouse CD84 (NP_038517.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q18PI6
Protein Weight	37kDa
Shipping	Ice bag