


Anti-Mouse Cd8A Monoclonal Antibody (Apc) For Flow Cytometry - P01731

Item Number: CM0004917



Introduction

Rabbit monoclonal antibody against mouse CD8a, directly conjugated with APC. Optimized for flow cytometry (FC) to study T-cell coreceptor function, CTL activity, and immune checkpoint research. Suitable for mouse immunology models.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD8a mAb
Remarks/Alias	Ly-2; Ly-B; Ly-35; Lyt-2
Species	Mouse
GeneID (Human)	925
GeneID	12525
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 28-196 of mouse CD8a (NP_001074579.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01731
Protein Weight	27kDa
Shipping	Ice bag