

Anti-Mouse Cd8A Monoclonal Antibody For Fc - P01731

Item Number: CM0004776



Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD8a (Lyt-2) for flow cytometry. Suitable for detecting cytotoxic T cells and NK cells in mouse samples in immunology and immuno-oncology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD8a mAb
Remarks/Alias	Ly-2; Ly-B; Ly-35; Lyt-2
Species	Mouse
GeneID (Human)	925
GeneID	12525
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-196 of mouse CD8a (NP_001074579.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01731
Protein Weight	27KDa
Shipping	Ice bag