


Anti-Mouse Cd8A Monoclonal Antibody For Fc - P01731

Item Number: CM0005302



Introduction

Rabbit anti-Mouse CD8a monoclonal antibody for flow cytometry. Recognizes the extracellular domain (aa 28-196) of mouse CD8a. Suitable for immune cell profiling and T-cell research. Calculated MW 27 kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD8a mAb
Remarks/Alias	Ly-2; Ly-B; Ly-35; Lyt-2
Species	Mouse
GeneID (Human)	925
GeneID	12525
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 28-196 of mouse CD8a (NP_001074579.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01731
Protein Weight	27kDa
Shipping	Ice bag