

Anti-Mouse Cd94 Monoclonal Antibody () For Flow Cytometry - O54707

Item Number: CM0005109

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD94, validated for flow cytometry. Essential tool for studying NK cell receptor function, immune self-tolerance, and tumor microenvironment interactions.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD94 mAb
Remarks/Alias	CD94
Species	Mouse
GenelD (Human)	3824
GenelD	16643
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 59-179 of mouse CD94 (NP_034784.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	O54707
Protein Weight	21kDa
Shipping	Ice bag