


Anti-Mouse Cd95/Fas Monoclonal Antibody For Fc - P25446

Item Number: CM0006057



Introduction

-conjugated rabbit anti-mouse CD95/FAS monoclonal antibody optimized for flow cytometry. Ideal for studying FAS receptor-mediated apoptosis and immune cell homeostasis in mouse models.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD95/FAS mAb
Remarks/Alias	lpr; APO1; APT1; CD95; TNFR6; Tnfrsf6
Species	Mouse
GeneID (Human)	355
GeneID	14102
Immunogen	Recombinant protein of mouse CD95/FAS
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P25446
Protein Weight	37KDa
Shipping	Ice bag