

Abflo 647 Anti-Mouse Cd98 (Slc3A2) Monoclonal Antibody For Flow Cytometry - P10852

Item Number: CM0005288



Introduction

ABflo 647-conjugated monoclonal antibody targeting mouse CD98 (SLC3A2) heavy chain. Validated for flow cytometry, enabling cell surface staining and amino acid transporter studies. Cross-reacts with mouse samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD98 mAb
Remarks/Alias	4F2; Cd98; Ly10; Mdu1; 4F2HC; Ly-10; NACAE; Ly-m10; Mgp-2hc
Species	Mouse
GeneID (Human)	6520
GeneID	17254
Immunogen	Recombinant protein corresponding to amino acids 105-526 of mouse CD98 (NP_032603.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P10852
Protein Weight	58kDa
Shipping	Ice bag