

Anti-Mouse Cx3Cr1 Monoclonal Antibody For Flow Cytometry - Q9Z0D9

Item Number: CM0012240

Introduction

Anti-Mouse CX3CR1 rabbit monoclonal antibody validated for flow cytometry. Recognizes mouse CX3C chemokine receptor 1 (fractalkine receptor). Suitable for immunology and inflammation research. SWISS: Q9Z0D9.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CX3CR1 Monoclonal Antibody (488 Conjugate)
Remarks/Alias	mCX3CR1
Species	Mouse
Gene ID (Human)	1524
Gene ID (Mouse)	13051
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of mouse CX3CR1 (NP_034117.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9Z0D9
Protein Weight	40 kDa
Shipping	Ice bag