

Anti-Mouse F4/80 Monoclonal Antibody (Biotin) For Flow Cytometry - Q61549

Item Number: CM0012772

Introduction

Biotin-conjugated rat anti-mouse F4/80 monoclonal antibody validated for flow cytometry. Specific for mouse F4/80 (Adgre1), a 102 kDa macrophage marker ideal for immune cell profiling.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rat anti-Mouse F4/80 mAb
Remarks/Alias	Emr1; Ly71; F4/80; Gpf480; TM7LN3; DD7A5-7; EGF-TM7
Species	Mouse
Gene ID (Mouse)	13733
Gene ID (Human)	2015
Immunogen	Recombinant protein of mouse F4/80.
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q61549
Protein Weight	102 kDa
Shipping	Ice bag