

Anti-Mouse G-Csf Rabbit Monoclonal Antibody For Wb, Elisa - P09920

Item Number: CM0007680

Introduction

Rabbit monoclonal antibody targeting mouse granulocyte colony-stimulating factor (G-CSF). Suitable for Western blot and ELISA. For research on hematopoiesis and granulocyte development. Reacts with mouse G-CSF.

[Learn More](#)

Parameter	Value
Product Name	Mouse G-CSF Rabbit mAb
Remarks/Alias	Csfg; G-CSF; MGI-IG; Mouse G-CSF
Species	Mouse
GeneID (Human)	1440
GeneID	12985
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-208 of mouse G-CSF (NP_034101.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P09920
Protein Weight	22kDa
Shipping	Ice bag