

Anti-Mouse Gsdma Rabbit Polyclonal Antibody - Gsdma

Item Number: CM0015854

Introduction

Rabbit polyclonal antibody targeting mouse Gsdma (Gasdermin A), raised against recombinant protein (aa 248-347). Validated for WB and ELISA; cross-reacts with rat. For cell death and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	Gsdma Rabbit pAb
Remarks/Alias	Gsdm, Gsdm1, H312E, Gsdma1, Gsdma
Species	Mouse
Gene ID (Human)	284110
Gene ID (Mouse)	57911
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 248-347 of mouse Gsdma (NP_067322.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
Protein Weight	49kDa
Shipping	Ice bag