


Anti-Mouse H2-Aa Rabbit Polyclonal Antibody For Wb And Elisa - P14434

Item Number: CM0023282



Introduction

Rabbit polyclonal antibody targeting mouse MHC class II alpha chain (H2-Aa), validated for WB and ELISA. Cross-reacts with human. Useful for immunological research on antigen presentation and adaptive immunity.

[Learn More](#)

Parameter	Value
Product Name	H2-Aa Rabbit pAb
Remarks/Alias	Ia1, H2Aa, Ia-1, H-2Aa, Aalpha, IAalpha, I-Aalpha, H2-Aa
Species	Mouse
GeneID	14960 (mouse)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-208 of mouse H2-Aa (NP_034508.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P14434
Protein Weight	28kDa
Shipping	Ice bag