

Anti-Mouse H2-Ab1 Rabbit Monoclonal Antibody For Wb - P06345

Item Number: CM0006308

Introduction

High-quality rabbit monoclonal antibody against mouse H2-Ab1, validated for WB, IF-P, IHC-P, and ELISA. Ideal for immunology research and MHC class II studies. 30kDa, cross-reacts with mouse.

[Learn More](#)

Parameter	Value
Product Name	H2-Ab1 Rabbit mAb
Remarks/Alias	IAb; Ia2; Ia-2; Abeta; H-2Ab; H2-Ab; Rmcs1; I-Abeta; H2-Ab1
Species	Mouse
GeneID (human)	474382
GeneID	14961
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 31-219 of mouse H2-Ab1 (NP_996988.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	P06345
Protein Weight	30kDa
Shipping	Ice bag