


Anti-Mouse Ifn-Beta Rabbit Monoclonal Antibody For Wb And Elisa - P01575

Item Number: CM0012605



Introduction

Rabbit monoclonal antibody against mouse interferon beta (IFN-beta), validated for Western blot and ELISA. Recognizes a secreted cytokine pivotal in antiviral defense, immune regulation, and neuronal protection. Suitable for immunology and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	IFN-beta Rabbit mAb
Remarks/Alias	Ifb; IFNB; Ifl1da1; IFN-beta
Species	Mouse
Gene ID (Human)	3456
Gene ID (Mouse)	15977
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-182 of mouse IFN-beta (NP_034640.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P01575
Protein Weight	22kDa
Shipping	Ice bag