

Anti-Mouse IgM Monoclonal Antibody For Wb, Ihc-P, Fc, Elisa - P01872

Item Number: CM0009351



Introduction

Recombinant rabbit monoclonal antibody targeting mouse immunoglobulin heavy constant mu (P01872) for Western blot, IHC-paraffin, flow cytometry, and ELISA. Validated in mouse samples, recognizes natural and immune IgM. Suitable for humoral immunity and B cell research.

[Learn More](#)

Parameter	Value
Product Name	Mouse IgM Rabbit mAb
Remarks/Alias	Igm; muH; Igh6; Igh-6; Igh-M
Species	Mouse
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 105-454 of mouse Ighm (P01872).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, FC, ELISA
Cross Reactivity	Mouse
SWISS	P01872
Protein Weight	50kDa/53kD
Shipping	Ice bag