

Anti-Mouse IgM Rabbit Monoclonal Antibody For Wb, If, Ihc, Elisa - P01872

Item Number: CM0005416



Introduction

Rabbit monoclonal antibody against mouse IgM constant region (Fc). Validated for WB, IF, IHC, ELISA. Cross-reacts with rat. Suitable for B cell research, BCR signaling, and humoral immunity studies.

[Learn More](#)

Parameter	Value
Product Name	Mouse IgM Rabbit mAb
Species	Mouse
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 105-454 of mouse IgM (P01872)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P01872
Protein Weight	50kDa
Shipping	Ice bag