


Anti-Mouse IL-12 P40 Monoclonal Antibody () For Intracellular Flow Cytometry - P43432

Item Number: CM0005254



Introduction

Rabbit monoclonal antibody against mouse IL-12 p40, conjugated to , validated for intracellular flow cytometry. Ideal for studying IL-12 in T cell, NK cell, and inflammatory responses.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse IL-12 p40 mAb
Remarks/Alias	p40; IL-12b; IL12p40; IL-12p40
Species	Mouse
Gene ID (Human)	3593
Gene ID	16160
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-335 of mouse IL-12 p40 (NP_001290173.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	P43432
Protein Weight	38kDa
Shipping	Ice bag