


Anti-Mouse IL-12 P40 Monoclonal Antibody For Intracellular Flow Cytometry - P43432

Item Number: CM0005194



Introduction

conjugated rabbit monoclonal antibody targeting mouse IL-12 p40, validated for intracellular flow cytometry. Recombinant immunogen ensures specificity for immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse IL-12 p40 mAb
Remarks/Alias	p40; IL-12b; IL12p40; IL-12p40
Species	Mouse
GeneID (Human)	3593
GeneID	16160
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-335 of mouse IL-12 p40(NP_001290173.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	P43432
Protein Weight	38kDa
Shipping	Ice bag