

Anti-Mouse IL1R1/Cd121A Monoclonal Antibody () For Flow Cytometry - P13504

Item Number: CM0005335

Introduction

-conjugated rabbit anti-mouse IL1R1/CD121a mAb, validated for flow cytometry. Targets mouse IL-1 receptor type 1 (67kDa, SWISS P13504); suitable for immunology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse IL1R1/CD121a mAb
Remarks/Alias	IL-iR; CD121a; CD121b; IL-1R1; IL1r-1; IL-1R-1; IL-1RT1; IL-1RT-1; IL-1R-alpha
Species	Mouse
GenelD (human)	3554
GenelD	16177
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P13504
Protein Weight	67kDa
Shipping	Ice bag