

Anti-Mouse IL4 Monoclonal Antibody (Apc) For Fc - P07750

Item Number: CM0005762

Introduction

APC-conjugated rabbit monoclonal antibody against mouse IL4 for intracellular flow cytometry. Specifically detects mouse IL4, a Th2 cytokine regulating B-cell activation, IgE class switching, and T-cell development.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse IL4 mAb
Remarks/Alias	IL-4; BSF-1
Species	Mouse
Gene ID (Human)	3565
Gene ID	16189
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-140 of mouse IL4 (NP_067258.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	P07750
Protein Weight	16 kDa
Shipping	Ice bag