

## Anti-Mouse Integrin Alpha 1/Cd49A Monoclonal Antibody () For Flow Cytometry - Q3V3R4

Item Number: CM0004661

### Introduction

-conjugated rabbit monoclonal antibody against mouse Integrin alpha 1/CD49a for flow cytometry applications. Validated in mouse samples. Predicted molecular weight 130 kDa.

[Learn More](#)

| Parameter        | Value                                                                                                                                   |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Mouse Integrin alpha 1/CD49a mAb                                                                                            |
| Remarks/Alias    | Vla1; CD49A; E130012M19Rik                                                                                                              |
| Species          | Mouse                                                                                                                                   |
| GeneID (Human)   | 3672                                                                                                                                    |
| GeneID           | 109700                                                                                                                                  |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 29-1141 of mouse Integrin alpha 1/CD49a (NP_001028400.2). |
| Source           | Rabbit                                                                                                                                  |
| Category         | Monoclonal Antibodies                                                                                                                   |
| Application      | FC                                                                                                                                      |
| Cross Reactivity | Mouse                                                                                                                                   |
| SWISS            | Q3V3R4                                                                                                                                  |
| Protein Weight   | 130kDa                                                                                                                                  |
| Shipping         | Ice bag                                                                                                                                 |