

Anti-Mouse Ki67 Monoclonal Antibody () For Flow Cytometry (Intracellular) - E9Pvx6

Item Number: CM0009789

Introduction

-conjugated rabbit monoclonal antibody targeting mouse Ki67, a proliferation marker. Validated for intracellular flow cytometry in mouse samples. Ideal for cell cycle and proliferation research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse Ki67 mAb
Remarks/Alias	Ki67; Ki-67; D630048A14Rik
Species	Mouse
GeneID (Human)	4288
GeneID	17345
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1816-2163 of mouse Ki67 (NP_001074586.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	E9PVX6
Protein Weight	351kDa
Shipping	Ice bag