

Anti-Mouse Ly-6A/E (Sca-1) Monoclonal Antibody For Fc - P05533

Item Number: CM0008116

Introduction

Rabbit anti-mouse Ly-6A/E (Sca-1) monoclonal antibody () validated for flow cytometry. Targets a GPI-anchored protein involved in T-cell activation. Suitable for mouse immunology research and stem cell studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse Ly-6A/E (Sca-1) mAb
Remarks/Alias	Ly-6A.2, Ly-6A/E, Ly-6E.1, Sca-1, Sca1, TAP
Species	Mouse
GeneID	110454
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-111 of mouse Ly-6A/E (Sca-1)(NP_034868.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P05533
Protein Weight	114kDa
Shipping	Ice bag