

Anti-Mouse Ly-6G/Ly-6C (Gr-1) Monoclonal Antibody () For Fc - P35461 / P0Cw02

Item Number: CM0009203

Introduction

-conjugated rat monoclonal antibody targeting mouse Ly-6G/Ly-6C (Gr-1) for flow cytometry. Enables reliable identification of granulocytes, monocytes, and myeloid-derived suppressors in immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rat anti-Mouse Ly-6G/Ly-6C (Gr-1) mAb
Remarks/Alias	Gr1; Gr-1; Ly-6G
Species	Mouse
Gene ID (Human)	8600
Gene ID	546644/17067
Immunogen	Recombinant protein of mouse Ly-6G/Ly-6C (Gr-1).
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P35461/P0CW02
Protein Weight	14kDa
Shipping	Ice bag