


Anti-Mouse Mhc Class Ii (I-A/I-E) Monoclonal Antibody For Flow Cytometry - P18468 / P06342 / Q3U060 / P14435

Item Number: CM0010241



Introduction

-conjugated rat anti-mouse MHC Class II (I-A/I-E) monoclonal antibody, validated for flow cytometry. Ideal for immune cell phenotyping, histocompatibility analysis, and diagnostic R&D.

[Learn More](#)

Parameter	Value
Product Name	Rat anti-Mouse MHC Class II (I-A/I-E) mAb
Remarks/Alias	Eb; Ia4; H2Eb; Ia-4; H-2Eb;IAb; Ia2; Ia-2; Abeta; H-2Ab; H2-Ab; Rmcs1; I-Abeta; Ia5; Ia-5; H-2Eb2; Iebeta2; A130038H09Rik; Ia1; H2Aa; Ia-1; H-2Aa; Aalpha; IAalpha; I-Aalpha
Species	Mouse
Gene ID	14969 / 14961 / 381091 / 14960
Immunogen	Recombinant protein of mouse MHC Class II (I-A/I-E)
Source	Rat
Category	Monoclonal Antibodies
Application	Flow Cytometry (FC)
Cross Reactivity	Mouse
SWISS	P18468 / P06342 / Q3U060 / P14435
Protein Weight	26kDa / 30kDa / 32kDa
Shipping	Ice bag