

Anti-Mouse Pd-L2/Cd273/B7-Dc Monoclonal Antibody () For Flow Cytometry - Q9Wul5

Item Number: CM0009971



Introduction

Rabbit monoclonal antibody against mouse PD-L2 (CD273/B7-DC), conjugated to for flow cytometry. Recognizes an immune checkpoint ligand that inhibits T-cell proliferation via PD-1 binding. Ideal for murine immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse PD-L2/CD273/B7-DC mAb
Remarks/Alias	Btdc; B7-DC; PD-L2; F730015O22Rik
Species	Mouse
Gene ID (Human)	80380
Gene ID	58205
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9WUL5
Protein Weight	28kDa
Shipping	Ice bag