

Anti-Mouse Placenta Growth Factor (Plgf) Polyclonal Antibody For Wb And Elisa - P49764

Item Number: CM0018882

Introduction

Rabbit polyclonal antibody targeting mouse Placenta Growth Factor (PIGF). Validated for WB and ELISA, suitable for angiogenesis and endothelial cell growth studies. Reacts specifically with mouse PIGF.

[Learn More](#)

Parameter	Value
Product Name	Mouse PIGF Rabbit pAb
Remarks/Alias	PIGF; Plgf; Mouse PIGF
Species	Mouse
Gene ID	18654
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-158 of mouse PIGF (NP_032853.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P49764
Protein Weight	17kDa
Shipping	Ice bag