

Anti-Mouse Podoplanin Monoclonal Antibody For Fc - Q62011

Item Number: CM0011113



Introduction

PE-conjugated rabbit monoclonal antibody against mouse Podoplanin (Gp38/T1alpha), validated for flow cytometry. Specifically detects mouse Podoplanin, useful for lymphatic and cancer research. SWISS: Q62011.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse Podoplanin mAb
Remarks/Alias	E11; T1a; Gp38; OTS-8; T1alpha; RANDAM-2; T1-alpha
Species	Mouse
GeneID (Human)	10630
GeneID	14726
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 21-133 of mouse Podoplanin (NP_034359.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q62011
Protein Weight	18kDa
Shipping	Ice bag