

Anti-Mouse Siglec5/Siglec-F Monoclonal Antibody () For Fc - Q920G3

Item Number: CM0009604

Introduction

-conjugated rabbit anti-mouse Siglec5 (Siglec-F) monoclonal antibody for flow cytometry. Detects sialic acid-binding lectin in mouse; suitable for identifying myeloid lineage cells in bone marrow and blood.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse Siglec5 mAb
Remarks/Alias	Siglec5; mSiglec-F
Species	Mouse
GeneID (Human)	8778
GeneID	233186
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 18-437 of mouse Siglec5 (NP_663556.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q920G3
Protein Weight	61 kDa
Shipping	Ice bag