

Anti-Mouse Tcr γ/δ Monoclonal Antibody - Tcr γ/δ

Item Number: CM0007809

Introduction

-conjugated Armenian Hamster anti-Mouse TCR γ/δ monoclonal antibody validated for flow cytometry. Detects mouse TCR γ/δ + T cells. Recombinant immunogen, specific for Mouse. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	anti-Mouse TCR γ/δ mAb
Remarks/Alias	Tcrdelta; TCRGV1S1; TCRGV2S1; TCRGV3S1; TCRGV5S3
Species	Mouse
Gene ID	110066/110067
Immunogen	Recombinant protein
Source	Armenian Hamster
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
Shipping	Ice bag