

Anti-Mouse Ythdc1 Polyclonal Antibody For Wb - E9Q5K9

Item Number: CM0024134

Introduction

Rabbit polyclonal antibody against mouse Ythdc1, validated in WB and ELISA. Targets amino acids 597-736. Suitable for m6A RNA binding and alternative splicing research.

[Learn More](#)

Parameter	Value
Product Name	Mouse Ythdc1 Rabbit pAb
Remarks/Alias	mKIAA1966; A730098D12Rik; Mouse Ythdc1
Species	Mouse
GeneID (Human)	91746
GeneID (Mouse)	231386
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 597-736 of mouse Ythdc1 (NP_808348.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	E9Q5K9
Protein Weight	85kDa
Shipping	Ice bag