

Anti-Mouse Zfp541 Rabbit Polyclonal Antibody - Zfp541

Item Number: CM0015222

Introduction

Rabbit polyclonal antibody targeting mouse Zfp541, suitable for WB and ELISA. Detects mouse and rat Zfp541. Immunogen: recombinant protein (aa 1-300). Predicted molecular weight: 141kDa/148kDa.

[Learn More](#)

Parameter	Value
Product Name	Zfp541 Rabbit pAb
Remarks / Alias	Znf541; EG666528; Zfp541
Species	Mouse
Gene ID (Human)	84215
Gene ID (Mouse)	666528
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of mouse Zfp541 (NP_001092747.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
Protein Weight	141kDa/148kDa
Shipping	Ice bag