

## Anti-Mouse/Rat Cd31/Pecam1 Monoclonal Antibody ( ) For Flow Cytometry - Q08481

Item Number: CM0009463

### Introduction

-conjugated rabbit monoclonal antibody against mouse/rat CD31/PECAM1. Validated for flow cytometry, ideal for leukocyte transendothelial migration, endothelial cell adhesion, and vascular biology research. Recognizes mouse and rat species.

[Learn More](#)

| Parameter        | Value                                                                                                                                    |
|------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Mouse/Rat CD31/PECAM1 mAb                                                                                                    |
| Remarks/Alias    | Cd31; Pecam; PECAM-1                                                                                                                     |
| Species          | Mouse                                                                                                                                    |
| GeneID (human)   | 5175                                                                                                                                     |
| GeneID           | 18613                                                                                                                                    |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 18-220 of mouse CD31 (NP_001292086.1). |
| Source           | Rabbit                                                                                                                                   |
| Category         | Monoclonal Antibodies                                                                                                                    |
| Application      | FC                                                                                                                                       |
| Cross Reactivity | Mouse, Rat                                                                                                                               |
| SWISS            | Q08481                                                                                                                                   |
| Protein Weight   | 69kDa/77kDa/80kDa/81kDa                                                                                                                  |
| Shipping         | Ice bag                                                                                                                                  |