

Anti-Mouse/Rat Cd90.1 Monoclonal Antibody For Flow Cytometry - P01831

Item Number: CM0010957

Introduction

-conjugated mouse monoclonal antibody specific for mouse/rat CD90.1 (Thy-1) membrane glycoprotein. Validated for flow cytometry (FC). Recognizes the GPI-anchored target involved in cell-cell interactions. Ideal for immunophenotyping and neuronal research.

[Learn More](#)

Parameter	Value
Product Name	Mouse anti-Mouse/Rat CD90.1 mAb
Remarks/Alias	T25; CD90; Thy-1; Thy1.1; Thy1.2; Thy-1.2
Species	Mouse
GeneID (Human)	7070
GeneID	21838
Immunogen	Recombinant protein of mouse/rat CD90.1
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse, Rat
SWISS	P01831
Shipping	Ice bag