

Anti-Mov10 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Hce1

Item Number: CM0012377

Introduction

Rabbit monoclonal antibody against MOV10 RNA helicase; validated for WB and ELISA; detects human, mouse, and rat MOV10; for studies on antiviral immunity, retrotransposition, and microRNA-mediated silencing.

[Learn More](#)

Parameter	Value
Product Name	MOV10 Rabbit mAb
Remarks/Alias	gb110; fSAP113; MOV10
Species	Human
GeneID (Human)	4343
GeneID	4343
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HCE1
Protein Weight	114kDa
Shipping	Ice bag