


Anti-Mpc2 Polyclonal Antibody For Wb, Ihc-P, Elisa - 095563

Item Number: CM0020155



Introduction

MPC2 rabbit polyclonal antibody (alias BRP44, SLC54A2) for WB, IHC-P, ELISA. Mitochondrial pyruvate carrier 2 mediates pyruvate uptake, balancing glycolysis and oxidative phosphorylation. Cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	MPC2 Rabbit pAb
Remarks/Alias	BRP44; SLC54A2; MPC2
Species	Human
GeneID (Human)	25874
GeneID	25874
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-60 of human MPC2 (NP_056230.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	095563
Protein Weight	14kDa
Shipping	Ice bag