


Anti-Mpdz Rabbit Polyclonal Antibody For Wb, Elisa - O75970

Item Number: CM0017927



Introduction

Rabbit polyclonal antibody targeting multiple PDZ domain protein (MPDZ/MUPP1), validated for WB and ELISA. Cross-reacts with human, mouse, rat. Ideal for synaptic plasticity and tight junction research.

[Learn More](#)

Parameter	Value
Product Name	MPDZ Rabbit pAb
Remarks/Alias	HYC2; MUPP1; MPDZ
Species	Human
GeneID	8777
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1100-1380 of human MPDZ (NP_003820.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75970
Protein Weight	222kDa
Shipping	Ice bag