


Anti-Mpg Rabbit Polyclonal Antibody (Ko Validated) For Wb, Elisa - P29372

Item Number: CM0024197



Introduction

Validated rabbit polyclonal antibody against MPG (AAG), a DNA-3-methyladenine glycosylase involved in base excision repair. Suitable for WB, ELISA with cross-reactivity to human, mouse, rat - P29372.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] MPG Rabbit pAb
Remarks/Alias	AAG; MDG; ADPG; APNG; Mid1; anpg; PIG11; PIG16; CRA36.1; PG
Species	Human
Gene ID	4350
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-293 of human MPG (NP_001015052.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29372
Protein Weight	33kDa
Shipping	Ice bag