


Anti-Mpp1 Polyclonal Antibody For Wb, Elisa - Q00013

Item Number: CM0028371



Introduction

High-quality rabbit polyclonal antibody targeting MPP1 (55 kDa erythrocyte membrane protein). Suitable for WB and ELISA, with cross-reactivity to human and mouse. Ideal for research on neutrophil polarity and membrane protein function.

[Learn More](#)

Parameter	Value
Product Name	MPP1 Rabbit pAb
Remarks/Alias	MRG1, PEMP, AAG12, EMP55, DXS552E, MPP1
Species	Human
Gene ID (Human)	4354
Gene ID	4354
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human MPP1 (NP_002427.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q00013
Protein Weight	52kDa
Shipping	Ice bag