

Anti-Mpp6 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nzw5

Item Number: CM0014335

Introduction

Rabbit polyclonal antibody targeting MPP6 (PALS2, VAM1) for WB and ELISA applications. Cross-reacts with human, mouse, rat. 61kDa membrane protein involved in scaffolding at cell junctions. Ideal for neuroscience and cancer research.

[Learn More](#)

Parameter	Value
Product Name	MPP6 Rabbit pAb
Remarks/Alias	MPP6; VAM1; p55T; VAM-1
Species	Human
GeneID (Human)	51678
GeneID	51678
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 220-420 of human MPP6 (NP_057531.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NZW5
Protein Weight	61kDa
Shipping	Ice bag