

Anti-Mpzl1 Polyclonal Antibody For Wb, Elisa - 095297

Item Number: CM0013095



Introduction

Rabbit polyclonal antibody against human MPZL1 (PZR), validated for Western blot and ELISA. Detects the 29 kDa cell surface receptor involved in ConA signaling and integrin-mediated motility, with cross-reactivity to mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	MPZL1 Rabbit pAb
Remarks/Alias	PZR; PZRα; PZRβ; PZR1b; MPZL1b; MPZL1
Species	Human
GeneID	9019
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 195-269 of human MPZL1 (NP_003944.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95297
Protein Weight	29kDa
Shipping	Ice bag