


Anti-Mras Rabbit Monoclonal Antibody For Wb, Elisa - O14807

Item Number: CM0003853



Introduction

Rabbit monoclonal antibody targeting human Mras (R-Ras3), raised against a synthetic peptide (aa 109-208). Validated for WB and ELISA, cross-reacts with human, mouse, and rat. Ideal for studying Ras-MAPK signaling in cancer and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	Mras Rabbit mAb
Remarks/Alias	NS11; M-RAS; RRAS3; R-RAS3; Mras
Species	Human
Gene ID	22808
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 109-208 of human Mras (O14807)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14807
Protein Weight	24kDa
Shipping	Ice bag