

Anti-Mrc2 Polyclonal Antibody For Wb, Elisa - Q9Ubg0

Item Number: CM0028544



Introduction

Rabbit polyclonal antibody against human MRC2 (CD280, UPARAP, ENDO180), validated for WB and ELISA. Detects endogenous MRC2 at ~167 kDa in human and mouse samples. Ideal for cancer and collagen remodeling research.

[Learn More](#)

Parameter	Value
Product Name	MRC2 Rabbit pAb
Remarks/Alias	CD280; UPARAP; CLEC13E; ENDO180; MRC2
Species	Human
GeneID (Human)	9902
GeneID	9902
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-300 of human MRC2 (NP_006030.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UBG0
Protein Weight	167kDa
Shipping	Ice bag