


Anti-Mre11 Polyclonal Antibody For Wb, Ip, Elisa - P49959

Item Number: CM0018165



Introduction

MRE11 Rabbit Polyclonal Antibody validated for WB, IP, ELISA. Detects endogenous MRE11 in Human and Rat samples. Targets double-strand break repair protein MRE11, a key component of the MRN complex essential for DNA damage response and genomic stability.

[Learn More](#)

Parameter	Value
Product Name	MRE11 Rabbit pAb
Remarks/Alias	ATLD; HNGS1; MRE11A; MRE11B; MRE11
Species	Human
GeneID (Human)	4361
GeneID	4361
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 290-485 of human MRE11 (NP_005582.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Rat
SWISS	P49959
Protein Weight	81kDa
Shipping	Ice bag