

Anti-Mre11 Rabbit Monoclonal Antibody For Wb, Elisa, Chip - P49959

Item Number: CM0005448

Introduction

Rabbit monoclonal antibody targeting human MRE11 (P49959), validated for WB, ELISA, and ChIP. Reacts with human and mouse. Critical for DNA double-strand break repair research.

[Learn More](#)

Parameter	Value
Product Name	MRE11 Rabbit mAb
Remarks/Alias	ATLD; HNGS1; MRE11A; MRE11B; MRE11
Species	Human
Gene ID (Human)	4361
Gene ID	4361
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 400-500 of human MRE11 (P49959).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA, ChIP
Cross Reactivity	Human, Mouse
SWISS	P49959
Protein Weight	81kDa
Shipping	Ice bag